



91%
overall
rating

Pac-Dent, Inc.
www.pac-dent.com

Description

Mighty Mini 90 is a single-use, disposable prophylaxis angle designed with a small, short head and slim profile to improve access and visibility during use. It is used with an ISO E-Type low-speed handpiece. This prophylaxis angle incorporates a proprietary beveled gear design for smooth running, reduced vibration and quiet operation. In addition, the ribbing on the external surface and the excellent flare of the cup make subgingival, lingual and interproximal access much easier and more efficient. It is latex-free. **Mighty Mini 90** is available in boxes of 144 prophylaxis angles with soft cups.

Indication

- For prophylaxis procedures

Unique Features

- Smooth-running, quiet beveled gear design
- Short head and slim profile for increased visibility and improved access
- Improved flare for better subgingival access
- Latex-free and ergonomic

“GOOD VISIBILITY WITH THE SMALL HEAD, SLIM PROFILE AND NARROW NECK.”

Consultants' Comments

- “A great prophylaxis angle to get to hard-to-reach areas.”
- “Great product - runs smoothly without any vibration, and good power level.”
- “The cups flared nicely and went subgingivally with ease - good interproximal access.”
- “There was minimal splatter.”
- “Stable without vibration - the prophylaxis angle stays on the motor.”


Evaluation Highlights

Mighty Mini 90 was evaluated by 17 consultants, with a total of 1,197 uses.


- Smooth-running, quiet beveled gear design.
- The smaller powerful head for increased visibility and improved access.
- Improved flare for better subgingival access.
- Ergonomic shape reduces hand and arm micro-traumas.
- 100% latex-free cup minimizes splatter.



Key Features:



Compared to Competitive Products:



Percentage of Consultants Who Would:

